

Pharmacological evaluation of bispecific antibodies: A novel method to quantify their *in vitro* activity in hematological samplesDaniel Primo<sup>1</sup>; Pilar Hernández<sup>1</sup>; Julián Gorrochategui<sup>1</sup>; María Luisa Vicente<sup>1</sup>; Cristina Gómez<sup>1</sup>; David Martínez-Cuadrón<sup>2</sup>; Pau Montesinos<sup>2</sup>; Joaquín Martínez-López<sup>3</sup>; Joan Ballesteros<sup>1</sup>.<sup>1</sup>Vivia Biotech, Madrid; <sup>2</sup>Hematology, Hospital Universitari i Politècnic La Fe de Valencia, Valencia, Spain; <sup>3</sup>Hospital Universitario 12 de Octubre, Madrid, Spain

## ABSTRACT

**Background and Aim:** The PharmaFlow automated flow platform has achieved 85% clinical correlation with AML samples with its novel Native Environment assay. Recently, novel Bi-specific antibodies (BsAbs) or analogous constructions that simultaneously bind a surface target on tumor cells and an associate TCR chain have been used as immunotherapy leading to T-cell activation and serial lysis of tumor cells that carry the target. The aim of the present study is develop and *in vitro* assay that incorporate multiples variables to better quantify the activity of bispecific antibodies and capture the interpatient variability.

**Material & Methods:** For this purpose, different fresh whole Bone Marrow (BM) or Peripheral Blood (PB) were tested with their corresponding BsAbs at 8 different concentrations in different time points (24h-144h). In this sense, we tested 31 AML BM samples (5 paired BM and PB) with the CD123xCD3 (Creative Biolabs) and 7 PB CLL and 3 B-ALL samples (one of them paired BM and PB) with Blinatumomab. When appropriate, basal quantification of TAA was performed by flow cytometry (FCM). The PharmaFlow platform efficiently count by FCM how many tumor cells are killed by every activated T-cells, here called effective E:T ratio. For each sample, 8-colour FCM staining was performed to simultaneously analyze the leukemic population, activated CD4 and CD8 T-cells and the residual normal cells.  $EC_{50}$  or  $E_{max}$  was calculated to evaluate potency or efficacy. Kinetics of activity was measured repeating the dose response curves in 3 different days.

**Results:** Most of the samples present both T-cell activation (CD25+) and an effective lysis of tumor cells after BsAbs exposure in a time and dose dependent manner, even starting with low basal E:T ratios (<1:100). For AML, basal quantification of CD123 by FCM density does not reflect a correlation with the *in vitro* response. By contrast, differences in T-cell cytotoxicity or leukemic immunoresistance were observed between samples in terms of  $EC_{50}$  or  $E_{max}$ , even more marked between CLL samples. The integration of effective E:T ratios,  $EC_{50}$ ,  $E_{max}$ , and kinetics allow us to generate an *in vitro* response model and select those samples with higher T-cell cytotoxicity after the different BsAbs exposure. Interestingly, many of the samples for all the BsAbs leave a significant proportion of live cells, even at the higher BsAb concentrations or with a remarkable expansion of activated T-cells that suggest the use of immunecheckpoint to unblock this immunoresistant status.

**Conclusion:** We have developed an automated flow cytometry assay for bispecific antibodies screening that keep intact both basal effector to target (E:T) ratios and Native environment using whole blood or bone marrow samples. In this context, the PharmaFlow platform selects different *in vitro* T-cytotoxicity effects across patients identifying best patient candidates for adoptive antitumor immunotherapy with BsAbs. The integration of Effective E:T ratios and pharmacological parameters better predict the *in vitro* response of BsAbs. Because of the high capacity of the PharmaFlow platform, additional antibodies constructions alone or in combinations with immunomodulatory agents could be tested to identify the better agents or immunotherapeutics combinations in hematological diseases.

## METHODS




Figure 1. Screening set-up and Workflow

## RESULTS

BsAb decrease leukemic cells and increase activated T-cells in a time and a concentration manner




Figure 2. Dose response curves to assess the CD3-CD123 bispecific antibody activity at different time points (72-96-120h) in AML samples. Upper panel displays leukemic cell depletion curves. The survival index (y-axis) ranges from 100% to 0% displaying the leukemic cell depletion after exposure to dose response CD3-CD123 bispecific antibody concentrations (x-axis). Bottom panel shows the simultaneous T-cell activation and proliferation along different time incubations. Absolute cell count of activated T-cells (y-axis) after CD3-CD123 bispecific antibody dose response concentrations.

## RESULTS


Simple Version Immune-Tumor Response  
How Activated (CD25+) T Cells Lead to Tumor Depletion??

Figure 3

## Activated T cells are the real drug: Effective E:T Ratios




Figure 4


New Method Effective E:T Ratios Captures I-O Activity  
Bispecific Abs

Figure 5: Representative example of the Effective T:E (y-axis) ratio from 8 AML samples (x-axis). This method clearly stratify patients with high (left side) vs low T-cell killing activity (right side).

Quantitative Pharmacology for Bispecific Antibodies  
Activity In Patient Samples

- EC50 tumor depletion (same T Cell proliferation)**
  - Very similar across most samples
  - When very low predicts patient may respond at low doses
  - When very high predicts resistant patient
- Effective E:T Ratio equivalent standard EC50**
  - Can be validated measuring dose responses with FACS sorted activated T Cells
  - High Effective E:T Ratios predicted sensitive patients
- E<sub>max</sub>**
  - E<sub>max</sub> near 100% required for a sensitive patient
  - High % resistant tumor cells suggest combination e.g. PD1
- Kinetics of response**

The integration of all these parameters quantifies the BsAb activity selecting cases with higher possibility to BsAb response.

Figure 6

## CONCLUSION

- We report a novel proprietary *ex vivo* automated flow cytometry assay for I-O drugs keeping intact both basal effector to target (E:T) ratios and native environment using whole BM or PB.
- The PharmaFlow platform selects different *in vitro* T-cytotoxicity effects across patients identifying best patient candidates for adoptive antitumor immunotherapy with BsAbs with the integration of Effective E:T ratios and pharmacological parameters ( $EC_{50}$  &  $E_{max}$ ): quantitative pharmacology of BsAbs in patient samples.
- Our findings are consistent with a model in which BsAb can enrich highly cytotoxic clonal T-subsets with Tumor-Specific Antigen in some patients.
- This assay enable evaluate multiple combinations with immunomodulators (PD1, CTLA-4, TIM-3, LAG-3) or BsAbs candidates for hematological diseases.
- Clinical trials should not exclude patient for low expression of TargetxCD3.
- New design of multi-specific antibodies from our new MOA are empowered by our screening of hundreds constructs *ex vivo*.
- CDx opportunity may increase substantially the clinical outcomes (ISTs).


Effective E:T ratios on the  
AML samples

Figure 7: Effective E:T ratios on AML samples (N=5) after the CD3-CD123 exposure (72h). Upper panel correspond to the BM and bottom panel to the PB compartment. Red Line correspond to the Leukemic cells and green line to the activated (CD25+) T-cells (CD4+ or CD8+). BM T-cells are better killers in 3/5 samples (#1, #4 and #5) consistent with BM immunosuppressed Tumor Infiltrated Lymphocytes (TILs).


Effective E:T ratios on the  
CLL samples

Figure 8: Effective E:T ratios on CLL samples (N=7) after the CD3-CD19 exposure (120h). Red Line correspond to the leukemic cells and green line to the activated (CD25+) T-cells (CD4+ or CD8+). Effective E:T ratios, AUC values,  $EC_{50}$  and  $E_{max}$  reflect both T-cell efficacy and B-cell immunoresistance identifying patients likely to receive blinatumomab therapy.


Effective E:T ratios on the  
ALL samples

Figure 9: Effective E:T ratios on ALL samples (N=4) after the CD3-CD19 exposure (72h). Red Line correspond to the Leukemic cells and green line to the activated (CD25+) T-cells (CD4+ or CD8+). Samples #1 and #2 correspond to paired BM and PB patient sample, being BM T-cells (#1) better killers than PB (#2).

## PD-1 increases T-cell activity




Figure 10: Absolute cell counts in a CLL sample for activated T-cells (left panels) and Live tumor cells (right panel) incubated with Blinatumomab in presence (green lines) or absence (blue lines) of the immunomodulatory inhibitor checkpoint PD1. Results show increased numbers of both activated T-cells (CD4 and CD8) in presence of PD1, enhancing the overall tumor killing.