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
* 1^a line CYT-IDR → complete remission & autologous transplant
• 2^a line Relapse 9 months, treated with CYT-IDR-FLU → Resistant, identified by test
• 3^a line Vidaza (Azacitidine) → progression
• 4^a line Vivia Test CIT-MIT-ETO → Complete Response 0.9% blast BM

METHODS

ExviTech® Platform


Figure 1

Plate setup. Eight different concentrations of each drug or drug combination is run for the used treatment protocols. The max concentration used is listed.

	1	2	3	4	5	6	7	8	9	10	11	12
A												
CYTARABINE	150µM											
IDARUBICIN	1.5µM											
FLUDARABINE	127.5µM											
MITOXANTRONE	10µM											
CLOFARABINE	74.8µM											
DAUNORUBICIN	30µM											
DMSO 0.4%	100%											
B												
AMSACRINE	3µM											
ETOPOSIDE	22.5µM											
6-THIOGUANINE	49.5µM											
DMSO 0.4%	100%											

Whole sample vs. Isolated Leukocytes: A. Correlation pairs showing differences among EC50 values from the same samples tested either as isolated leukocytes or whole sample. Error bars show the CI's of the estimated parameter. B. Dose-response curves for IDA and Cyta for the selected samples in both conditions, showing similar results for Cytarabine but very different for Idarubicin.



Figure 2

Data Analysis: performed using the population approach using NONMEM 7.2:


- population PB Hill-based modelling of the ex vivo response vs concentration data in monotherapy (fig.2), 95% confidence interval of estimated parameters determined by bootstrapping over 1000 simulations.
- Surface interaction modelling and simulations to estimate the interaction parameter (α) as well as the corresponding confidence interval. α parameter is a measurement of synergism (>0), additivity (0) or antagonism (<0). Greco et al. 1995. Pharmacol Rev June 1995 47:331-385

RESULTS

Dose-responses from 377 patient samples to Cytarabine (A) and Idarubicin (B). The Survival Index (y-axis) ranges from 100% to 0 displaying the selective AML cell depletion. Median response shown in red. For CYT 40% patient samples have resistant cells left alive at 48 h. IDA eliminates all cells within this timeframe.

Interaction index

Distribution of the interaction index calculated for the samples in the study. The central tendency was towards an additive or weak synergistic behavior although clear synergy occurred in many cases.


Pharmacological Population Parameters

SINGLE DRUG ex vivo PHARMACOLOGY					
DRUG	N	Efficacy (E _{max}) % Survival	Potency (EC ₅₀) µM	IPV-E _{max} %	IPV-EC ₅₀ %
CYT	377	8.5 1.3	2.87 8.1	645 9.1	161 8.4
IDA	377	0 -	0.032 9.4	ne	207 8.4

Individual drug typical and relative standard error values. Inter-patient variability (IPV) expressed as CV(%); ne, not estimated

CORRELATION WITHOUT CYTOGENETICS

Figure 6 Logistic additive model of ex vivo CYT-IDA vs Clinical Outcome

- A generalized binary logistic additive model was used to explore nonparametric relationships between either the fitted pharmacologic parameters and processed response values and the dichotomized clinical response (resistant patient [PR or PD after induction] vs. sensitive patient [CR or CRi after induction]).
- Both linear dependence and nonlinear dependence structures were evaluated for available PD parameters (cytarabine E₀, EC₅₀, I_{max} and sigmoidicity, idarubicin E₀, EC₅₀ and sigmoidicity, and the interaction parameter (α)) as well as integrated terms given by the calculation of the area under the curve (AUC) for both cytarabine and idarubicin and the volume under the surface (VUS) from the interaction analysis from the combination.
- All linear terms were non-significant. Results using individual parameters were improved by the AUCs of the modelled effect-concentration curves of both, idarubicin and, particularly, Cytarabine which showed good predictive properties. In a lower magnitude, No significance though was observed for the interaction parameter.
- Using a criterion based on equalizing the predictive values (PV+ and PV-) to set the cut point which defines positive and negative test results is a reasonable approach to prioritize specificity over sensitivity in an objective and reproducible fashion.

Figure 7 Correlation results summary

Clinical outcome		Subtotal		Key clinical indicators: overall prediction 81.7% & NPV 81.9%
RESISTANT	SENSITIVE	26	18.3%	
21	5	26	18.3%	
14.8%	3.5%			
21	95	116	81.7%	
14.8%	66.9%			
Sensitivity %	Specificity %			
50.00	95.00	81.69		
42	100	N	142	
29.6%	70.4%		100.0%	

CONCLUSIONS

- Good predictive capabilities were identified for dose-effect area under the curve variables.
- No statistical significance with the clinical outcome was found for the interaction index from the drugs combination analysis.
- This novel test alone, without any clinical information, is able to predict the clinical response to Ida+Ara-C induction with 82% correlation with clinical outcome in 142 patient samples, significantly higher than the current clinical response rate of 70%. The test did not properly identify 21/142 that were clinically resistant and the test predicted as sensitive. This mismatched subgroup mimics the problems from molecular markers where a resistant clone present in a minority of leukemic cells cannot be detected yet drives the patient response. However, this group mismatch does not prevent a good correlation with the test predicted outcomes.
- This novel test including cytogenetic risk profiles for each patient is able to predict the clinical response to Ida+Ara-C induction with 90% correlation with clinical outcome. Cytogenetic information enables correct identification of resistant patients predicted sensitive by the ex vivo results (from 21/42 to 32/32 resistant patients).
- This PM test could improve substantially AML 1st line treatment outcomes.
- We are preparing an interventional clinical trial with PETHEMA.

ACKNOWLEDGEMENTS

Special Thanks to the Patients and Hospitals for Providing the Samples (listed alphabetically)

- Complejo Hospitalario de Jaén, JAEN
- Hospital Ramón y Cajal, MADRID
- Complejo Hospitalario Xeral Cies de Vigo, VIGO
- Hospital Universitario Central de Asturias, OVIEDO
- Hospital Carlos Haya, MÁLAGA
- Hospital Universitario de Canarias, TENERIFE
- Hospital Universitario San Cecilio, GRANADA
- Hospital Universitario General de Castellón, CASTELLÓN
- Hospital Universitario Gregorio Marañón, MADRID
- Hospital Universitari i Politècnic La Fe, VALENCIA
- Hospital Universitario Infanta Leonor, MADRID
- Hospital Universitario Lucas Augusti, LUGO
- Hospital Universitario Príncipe de Asturias, MADRID
- Hospital Universitario Reina Sofía, CÓRDOBA
- Hospital Universitario Virgen Macarena, SEVILLA