

THE PHARMACOLOGICAL PROFILES OF CYTOTOXIC AML TREATMENTS EX VIVO IDENTIFIES SENSITIVE VS RESISTANT LEUKEMIC CELLS

Joan Ballesteros Nobell¹, Pau Montesinos², Federico Moscardó², David Martínez Cuadrón², Miguel A. Sanz², Joaquín Martínez López², Jaime Pérez de Oteyza⁴, Raimundo García Boyero⁵, Josefina Serrano⁶, Pascual Fernández⁷, Pilar Herrera⁸, Ángeles Fernández⁹, Arancha Alonso¹⁰, Concepcion Bethancourt¹¹, Ataulfo Gonzalez¹², Pilar Hernández-Campo¹, Julian Gorrochategui¹, Belén Liébana¹, Iñaki F. Trocóniz¹³, Teresa A. Bennett¹

¹Vivia Biotech, Madrid, ²Haematology, Hospital Universitari i politècnic La Fe, Valencia, ³Haematology, Hospital Universitario Doce de Octubre, Madrid, ⁴Haematology, Hospital de Madrid Norte Sanchinarro, Madrid, ⁵Haematology, Hospital Universitario General de Castellón, Castellón, ⁶Haematology, Hospital Universitario Reina Sofía, Córdoba, ⁷Haematology, Hospital General Universitario de Alicante, Alicante, ⁸Haematology, Hospital Universitario Ramón y Cajal, ⁹Haematology, Complejo Hospitalario Xeral Cies de Vigo, Vigo, ¹⁰Haematology, Hospital Quirón Madrid, ¹¹Hospital Regional Universitario Carlos Haya, Málaga, ¹²Haematology, Hospital Clínico San Carlos, Madrid, ¹³Universidad de Navarra, Pamplona, Spain

ABSTRACT

Background: To aid in the identification of effective treatments for individual patients, ex vivo assays for detecting cell death inducible by drugs for hematological malignancies have been in development for over 20 years. We have developed a novel approach incorporating 4 key innovations: incubating drugs in whole bone marrow sample without isolating leukocytes, using flow cytometry enables identification of the malignant cells selectively, an automated flow cytometry-based platform (ExviTech) decreases errors and enables full pharmacological characterization, and analyzing the data using pharmacodynamic population models.

Aim: Derive the ex vivo pharmacological profiles across the AML patient population of single drugs and combination treatments as a tool for individualized treatment selection.

Patients and Methods: Bone-marrow samples from 160 patients diagnosed with AML were sent to Vivia from 24 hospitals in Spain within 24h. Plates incubated for 48-hours prior to analysis with ExviTech. Percentage of leukemic cell death was determined labeling with monoclonal antibodies and AnnexinV-FITC. Survival index is computed for each drug, the lower the survival index, the more effective the drug. Dose-response curves of cytarabine, idarubicin, daunorubicine, etoposide, mitoxantrone, fludarabine, clofarabine, and 6-thioguanine were measured in 160 samples. The added benefit of combining these drugs into 12 combination treatments was assessed by measuring their synergy in each individual patient. In 39 patients treated with CYT IDA we had clinical data of response, and then we performed a blinded interpretation of this in vitro test by an expert hematologist, to predict the clinical response based in this test result.

Results: There was a large range of interpatient variability in the response to a single drug and even larger in the synergism between drugs. Population Pharmacological Profiles for two individual patients are shown on the figure 6. Relative drug potency in terms of percentile ranking within the population is shown in the left panel from 0 (weakest) to 100 (most potent). Green lines show individual patient potency relative to the population ranking, with confidence intervals (CI). 3rd column lists when a drug leaves a significant % of leukemic cells alive, potential resistant clones. Synergism value for an individual patient in each combination is shown in green, with CI as parallel dotted green lines. Representation of the Pharmacological Profile of an individual patient sample quickly identifies extreme values, when a drug or combination is very sensitive (rightward shift green lines, green boxes) or very resistant (leftward shift green lines, red boxes). These representations lead to clear guidelines in >90% samples, and based on hematologist's interpretation of these guidelines show a clinical correlation with clinical responses to CYT-IDA of 84%.

METHODS

ExviTech® Platform





PLATE SETUP

Eight different concentrations of each drug or drug combination is run for the used treatment protocols. The max concentration used is listed

	1	2	3	4	5	6	7	8	9	10	11	12
A	CYTARABINE 150µM											
B	IDARUBICIN 1.5µM											
C	FLUDARABINE 127.5µM											
D	MITOXANTRONE 10µM											
E	CLOFARABINE 74.8µM											
F	PANOBIPOSTAT 30µM											
G	DAUNORUBICIN 3µM											
H	ETOPOSIDE 225µM											
I	6-THIOGUANINE 100µM											
J	CYCLOPHOSPHAMID 100µM											
K	MEPHALAN 15µM											
L	5-AZACYTIDINE 50µM											


Whole sample vs. Isolated Leukocytes: A. Dose-response curves for IDA and CYT in isolated leukocytes and whole sample. Data, from sample 6 below, displays a log difference in the EC50s for IDA, but equal results for CYT. B. The EC50 (y-axis) of the whole sample and the isolated leukocyte fraction from 9 patient samples with cytarabine. C. EC50 of the same samples to idarubicin.

Data Analysis: performed using the population approach using NONMEM 7.2.: population PD modelling of the ex vivo response vs concentration data in monotherapy (fig.1), establishing for each patient the 95% prediction intervals (PI) of the isobologram from each individual parameter (fig.4) computation of the combination index using raw data descriptors from combination experiments. Chou and Talalay, 2010. *Cancer Research* 70: 440-446.



RESULTS


Figure 1 Population PD modelling of the ex vivo response



A. Median curves from 8 drugs. The Survival Index (y-axis) ranges from 100% to 0 displaying the selective AML cell depletion calculated with PKPD Population Models. B. The curves from A. overlaid on 125 individual dose-response curves to fludarabine.


	IDA	DAU	MIT	ETO	CYT	FLU	CLO	THI
N	125	109	110	110	125	125	122	86
EC50 µM	0.1	0.6	0.2	14.6	2.3	1.4	0.9	62.2

Figure 2 Individual dose response curves per drug and their median DR curve



Dose-response analysis was completed for individual drug in 86-125 AML patient bone marrow samples. The Survival Index (y-axis) ranges from 100% to 0 displaying the selective AML cell depletion calculated with PKPD Population Models. The grey lines display each individual response with the median response shown in red.

Figure 4 Drug interaction description



Dose-response of 2 samples to Cyt alone (solid line) and Cyt+Ida (dashed line). A displays synergism; B an additive response. C. The Combination Index (CI): Synergistic CI<1, Additive CI=1 or Antagonistic CI>1.


Table 1 Pharmacological Population Parameters

SINGLE DRUG PHARMACOLOGY											
DRUG	N	Efficacy (E _{max}) % Survival	Potency (EC ₅₀) µM	IPV-E _{max}	IPV-EC ₅₀	Typical	RE	Typical	RE	Typical	RE
IDA	125	0*	-	0.106	0.016	ne	157	0.496	0.14		
DAU	109	0*	0	0.592	0.13	ne	135	0.17			
MIT	110	0.5	0.2	0.233	0.22	1.2	0.69	221	0.14		
ETO	110	0.1	0.1	18.5	0.13	0.3	0.41	141	0.17		
CYT	125	11.8	4	2.28	0.13	32	0.21	105	0.25		
FLU	125	15.0	3.2	1.43	0.19	25	0.24	113	0.41		
CLO	122	29.0	5.1	0.92	0.019	36	0.2	142	0.25		
THIO	86	0*	-	62.2	0.23	ne	204	0.27			

SYNERGY COMBINATION TREATMENT											
Combination	N	Median	SE	CI							
CYT+IDA	99	0.496	0.14								
CYT+FLU	82	0.548	0.24								
CYT+IDA+FLU	69	0.441	0.13								
CYT+DAU	65	0.743	0.26								
CYT+MIT	24	0.533	0.20								
CYT+DAU+FLU	8	0.677	0.27								
CYT+MIT+FLU	14	0.317	0.09								
CYT+IDA+ETO	26	0.535	0.16								
CYT+IDA+ETO	22	0.551	0.23								
CYT+MIT+ETO	13	0.421	0.32								
CYT+THIO	7	0.617	0.35								
CYT+CLO	53	0.552	0.44								


Individual drug typical and random error values (left). Inter-patient variability (IPV) expressed as CV(%); Synergism (right) using the CI. *, estimate not significantly different from 0; ne, not estimated

Figure 3 Pharmacologic Profile of ex vivo response to each drug.



Monotherapy profiles interpretation

Figure 5 Population synergy profiles



Synergy profiles of 12 drug combination treatments

ACKNOWLEDGEMENTS

Special Thanks to the Patients and Hospitals for Providing the Samples

- ✦ Complejo Hospitalario de Jaén, JAEN
- ✦ Hospital Universitario Central de Asturias, OVIEDO
- ✦ Hospital Carlos Haya, MÁLAGA
- ✦ Hospital Clínico San Carlos, MADRID
- ✦ Hospital de la Santa Creu i Sant Pau, BARCELONA
- ✦ Hospital de Madrid Norte Sanchinarro, MADRID
- ✦ Hospital Doce de Octubre, MADRID
- ✦ Hospital General Universitario de Alicante, ALICANTE
- ✦ Hospital General Universitario de Navarra, PAMPLONA
- ✦ Hospital General Universitario de Valencia, VALENCIA
- ✦ Hospital General Universitario de Zaragoza, ZARAGOZA
- ✦ Hospital General Universitario de Murcia, MURCIA
- ✦ Hospital General Universitario de Granada, GRANADA
- ✦ Hospital General Universitario de Sevilla, SEVILLA
- ✦ Hospital General Universitario de Valladolid, VALLADOLID
- ✦ Hospital General Universitario de Salamanca, SALAMANCA
- ✦ Hospital General Universitario de Burgos, BURGOS
- ✦ Hospital General Universitario de Cantabria, SANTANDER
- ✦ Hospital General Universitario de La Rioja, LOGROÑO
- ✦ Hospital General Universitario de Navarra, PAMPLONA
- ✦ Hospital General Universitario de Valencia, VALENCIA
- ✦ Hospital General Universitario de Zaragoza, ZARAGOZA
- ✦ Hospital General Universitario de Murcia, MURCIA
- ✦ Hospital General Universitario de Granada, GRANADA
- ✦ Hospital General Universitario de Sevilla, SEVILLA
- ✦ Hospital General Universitario de Valladolid, VALLADOLID
- ✦ Hospital General Universitario de Salamanca, SALAMANCA
- ✦ Hospital General Universitario de Burgos, BURGOS
- ✦ Hospital General Universitario de Cantabria, SANTANDER
- ✦ Hospital General Universitario de La Rioja, LOGROÑO

CONCLUSIONS

- By testing the drugs used in the treatment protocols for AML directly on patient samples, a pharmacological based model could be developed to infer drug resistance or sensitivity, patient by patient.
- Similarity, testing could be used as a companion diagnostic to identify subsets of patients for which specific cytotoxic drugs or targeted therapies would be effective.
- The Pharmacological Profiles could be used personalize treatment for individual patients.
- Correlation of this ex vivo sensitivity with the clinical efficacy is currently being performed in a study under the supervision of the PETHEMA group.